

An Automated Bioreactor Sampling Solution for Assuring On-line PAT Analytical Fidelity


2014 IFPAC Annual Meeting

*William Miller, Lucas Schimmelpfenning,
Paul Strand & Michael Biksacky*

Outline

- **Automated Bioreactor Sampling**
- **Seg-Flow Technology Overview**
- **Analytical Performance Case Studies**
 - Methodology & Acceptance Criteria
 - YSI 2700 Biochemistry Analyzer
 - Nova BioProfile 400 Analyzer
 - Vi-CELL XR Cell Analyzer
- **Summary**

Flownamics, Inc.



- **Global PAT solutions provider for bio & other processing industries**
 - Upstream process focus
- **Enabling technologies for on-line process analytics & data management**
- **Headquartered in Madison, WI USA**

Automated Bioreactor Sampling

- **Criteria**


- not compromise the bioreactor's sterile environment
- establish connectivity of process systems
- facilitate a scale-independent strategy
- seamlessly integrate the real-time data into the process management system
- **provide rapid and precise analysis that performs as good as or better than the manual off-line analytical method**

Technology Overview

SEGFLOW®



Seg-Flow Technology



Case Studies

Methodology

- **Objectives:**

- Does the Seg-Flow-integrated instrument perform as good as or better than the manufacturer's precision specifications?
 - Is the analytical fidelity preserved with the Seg-Flow system?
- Is the Seg-Flow automated on-line analytical method comparable to the manual off-line analytical method?

Methodology

- **Evaluation:**

- Evaluate the integrated analytical performance of three analyzers commonly used with bioreactor culture monitoring.
 - YSI® 2700 Select Biochemistry Analyzer
 - Nova® BioProfile® 400 Analyzer
 - Vi-Cell® XR Cell Analyzer

Methodology

- **General Scheme:**
 - Precision Evaluation
 - Within run evaluation per manufacturer's specifications
 - Comparability Evaluation
 - Evaluate $\geq 50\%$ of the instrument's measurement range
 - CDM/reagent standards used in lieu of live culture
 - Assure QC of analyte concentrations
 - 0.25 – 5.0L WV: serial dilutions to attain measurement ranges

Methodology

- **General Scheme:**

- Analytical instruments QC'd prior to evaluation
 - Manufacturer's linearity, QC standards used
- Seg-Flow/single instrument integration
 - Sample cycle – purge, analysis & system cleaning
 - precision & comparability studies
 - Manual sample analysis performed ≤ 5 minutes of Seg-Flow system analysis
 - comparability study

Acceptance Criteria

- **Performance standards based on:**
 - Instrument manufacturer's precision specifications
 - Accepted practices and standards

Acceptance Criteria

- **Precision:**
 - 2-point linearity check
 - Coefficient of variation (%CV) (YSI/Nova)
 - **% CV ≤ the manufacturer's within run specification**
 - compares the dispersion or variation in groups of measurements²
 - $\delta/\mu \times 100\%$
 - Concentration average accuracy (Vi-CELL XR)
 - **Average accuracy within $\pm 3.0\%$ of reference standard**

Acceptance Criteria

- **Accuracy:**
 - Qualitative evaluation of Seg-Flow System
 - **No analytical errors due to Seg-Flow sample delivery**
 - Ensure prescribed sample volume and timing are achieved

Acceptance Criteria

- **Comparability:**

- **Linear Regression Analysis**

- Determine statistical relationship of two analytical methods
- Correlation coefficient
 - **$R \geq 0.98^3$**
 - Strong positive linear correlation should exist between the Seg-Flow (automated) & manual analytical methods
- Slope
 - **95% CI should include the value of 1.0**
 - **slope = 1.0 (perfect)**
- Intercept
 - **95% CI should include the value of 0.0**
 - **intercept = 0 (perfect)**

Seg-Flow Integration: YSI® 2700 Biochemistry Analyzer



Analytes

- D-Glucose
- L-Lactate
- L-Glutamine
- L-Glutamate

Seg-Flow Integration: YSI® 2700 Biochemistry Analyzer





YSI 2700 Measurement Ranges & Precision Specifications				
Analyte	Measurement Range	Evaluated Analytical Range	CV (%)	Sample size (n)
D-Glucose	0 – 25.0 g/L	0.5 – 15.0 g/L	2.0	10
L-Lactate	0 - 2.7 g/L	0.2 – 5.0 g/L	2.0	10

Seg-Flow/YSI 2700 Results

Seg-Flow/YSI 2700 within run precision					
Analyte	Number of Samples Analyzed	Theoretical concentration (g/L)	Measured concentration (g/L) ($\mu \pm \delta$)	CV (%)	YSI CV Spec. (%)
D-Glucose	10	2.00	2.03 \pm 0.04	1.98	2.00
	10	10.00	10.01 \pm 0.20	1.96	2.00
L-Lactate	10	0.70	0.64 \pm 0.01	1.56	2.00
	10	2.70	2.74 \pm 0.05	1.94	2.00

- **YSI 2700 precision linearity demonstrated**
 - %CV acceptance criteria met
- **Seg-Flow sample delivery accuracy achieved**
 - No analytical errors due to sample delivery

Seg-Flow/YSI 2700 Results



Seg-Flow/YSI 2700 Results

Seg-Flow/YSI 2700 Statistical Comparability

Analyte	FISP Sampling Probe			Dip Tube		
	R	Slope (95% CI)	Intercept (95% CI)	R	Slope (95% CI)	Intercept (95% CI)
D-Glucose	1.00	1.007 (0.994 to 1.020)	0.0048 (-0.1087 to 0.1183)	1.00	0.999 (0.983 to 1.014)	-0.0385 (-0.1710 to 0.0940)
L-Lactate	1.00	0.998 (0.985 to 1.010)	0.028 (-0.007 to 0.062)	1.00	0.992 (0.977 to 1.007)	0.013 (-0.028 to 0.053)
L-Glutamate	1.00	0.985 (0.973 to 0.997)	0.0148 (-0.0150 to 0.0445)	1.00	0.983 (0.967 to 1.000)	-0.0385 (-0.1710 to 0.0940)
L-Glutamine	1.00	0.988 (0.969 to 1.007)	0.0411 (-0.0161 to 0.0984)	1.00	0.988 (0.969 to 1.007)	0.0411 (-0.0161 to 0.0984)

- **Comparability demonstrated for Seg-Flow & manual analytical methods**
 - Acceptance criteria met:
 - $R \geq 0.98$
 - Slope & Intercept within 95% CI
 - Irrespective of sampling mechanism used

Seg-Flow Integration: Nova BioProfile 400



Analytes

- Glucose
- Lactate
- Glutamine
- Glutamate
- Ammonium
- pO₂
- pCO₂
- pH
- Potassium
- Sodium
- Osmolality

Seg-Flow Integration: Nova BioProfile 400

Nova BP 400 Measurement Ranges & Within Run Precision Specifications





Analyte	Measurement Range	Evaluated Analytical Range	CV (%)	Sample size (n)
D-Glucose	0.2 – 15.0 g/L	0.5 – 15.0 g/L	5.0	20
L-Lactate	0.2 – 5.0 g/L	0.2 – 5.0 g/L	5.0	20
L-Glutamate	0.2 – 6.0 mmol/L	0.2 – 5.0 mmol/L	5.0	20
pO ₂	0 – 800 mmHg	170 – 230 mmHg	5.0	20
pCO ₂	3 – 300 mmHg	18 – 50 mmHg	5.0	20

Seg-Flow/Nova 400 Results




Seg-Flow/Nova 400 within run precision					
Analyte	Number of Samples Analyzed	Theoretical concentration	Measured concentration ($\mu \pm \delta$)	CV (%)	Nova CV Spec. (%)
D-Glucose	20	8.0 g/L	8.1 \pm 0.1	1.6	5.0
	20	15.0 g/L	15.0 \pm 0.4	2.5	5.0
L-Lactate	20	2.5 g/L	2.4 \pm 0.1	2.1	5.0
	20	5.0 g/L	4.7 \pm 0.1	2.9	5.0
L-Glutamate	20	2.5 mmol/L	2.3 \pm 0.1	2.7	5.0
	20	4.5 mmol/L	4.3 \pm 0.1	2.9	5.0
pO ₂	10	185 mmHg	185.7 \pm 2.7	0.9	5.0
	10	220 mmHg	220.5 \pm 4.5	2.0	5.0
pCO ₂	10	18 mmHg	18.2 \pm 0.3	1.9	5.0
	10	48 mmHg	47.5 \pm 0.8	1.8	5.0

- **Nova 400 precision acceptance criteria met**
 - Precision linearity demonstrated
- **Seg-Flow sample delivery accuracy achieved**
 - No analytical errors due to sample delivery

Seg-Flow/Nova 400 Results



Seg-Flow/Nova 400 Results



Seg-Flow/Nova 400 Results

Seg-Flow/Nova BP 400 Statistical Comparability

Analyte	FISP Sampling Probe			Dip Tube		
	R	Slope (95% CI)	Intercept (95% CI)	R	Slope (95% CI)	Intercept (95% CI)
D-Glucose	0.99	0.994 (0.928 to 1.059)	0.260 (-0.346 to 0.867)	0.99	1.024 (0.966 to 1.082)	0.027 (-0.508 to 0.563)
L-Lactate	0.99	0.942 (0.874 to 1.009)	-0.126 (-0.337 to 0.086)	0.99	0.998 (0.897 to 1.099)	-0.332 (-0.680 to 0.016)
L-Glutamate	0.99	0.994 (0.856 to 1.131)	-0.003 (-0.383 to 0.376)	0.99	0.986 (0.868 to 1.103)	-0.007 (-0.341 to 0.327)
L-Glutamine	1.00	1.048 (0.984 to 1.113)	-0.106 (-0.338 to 0.126)	0.99	1.053 (0.951 to 1.156)	-0.218 (-0.536 to 0.099)
Ammonium	1.00	1.044 (0.994 to 1.093)	-0.003 (-0.171 to 0.165)	1.00	1.028 (0.993 to 1.064)	0.058 (-0.052 to 0.168)
pO2	-	-	-	1.00	0.924 (0.816 to 1.032)	13.98 (-6.84 to 34.80)
pCO2	-	-	-	0.98	1.066 (0.842 to 1.290)	-3.97 (-13.50 to 5.57)

- **Comparability demonstrated for Seg-Flow & manual analytical methods**
 - Acceptance criteria met:
 - $R \geq 0.98$
 - Slope & Intercept within 95% CI
 - Irrespective of sampling mechanism used

Seg-Flow Integration: Vi-CELL® XR Cell Analyzer




Analytes

- VCC
- TCC
- % Viability
- Total cell count
- Viable cell count
- μ Cell Diameter
- μ Compactness
- Aggregation Rate
- Cell Imaging

- 12 hour test duration w/ 30 minute sample frequency
- # samples represent typical 2 - 4 week cell culture sampling (1 - 2/day)
- Cell concentration calibration beads used for analysis

Seg-Flow/Vi-CELL XR Results

Mean TCC, VCC and % Viability



- **Vi-CELL XR precision acceptance criteria met**
 - Concentration Average Accuracy: $\pm 3.0\%$ (n = 20)
 - $\pm 2.4\%$ average concentration accuracy for VCC (n=25)
 - $\pm 2.0\%$ average concentration accuracy for TCC (n = 25)
 - $\pm 0.5\%$ difference observed for % viability (n = 25)
- **Seg-Flow sample delivery accuracy achieved**
 - No analytical errors due to sample delivery

Seg-Flow/Vi-CELL XR Results

Seg-Flow/Vi-CELL XR Statistical Comparability			
Analyte	R	Slope (95% CI)	Intercept (95% CI)
VCC	0.98	1.047 (0.982 to 1.112)	-0.01806 (-0.05543 to 0.01931)
TCC	0.98	1.034 (0.982 to 1.085)	-0.01040 (-0.04090 to 0.02009)

- **Comparability demonstrated for Seg-Flow & manual analytical methods**
 - Acceptance criteria met:
 - $R \geq 0.98$
 - Slope & Intercept within 95% CI

Summary

- **Analytical fidelity (precision) preserved for each Seg-Flow-integrated analyzer**
 - YSI 2700 Biochemistry Analyzer
 - Nova BioProfile 400 Analyzer
 - Vi-CELL XR Cell Analyzer
- **Seg-Flow automated and manual analytical methods are statistically comparable**
 - Fully automated system for delivering precise, reliable analyses
 - cell parameters, nutrients, metabolites and product
- **Enabling on-line PAT solution for real-time bioreactor culture monitoring**
 - Achieve deeper process understanding & increase process efficiency

Acknowledgements



Jayson Preston
Ashley Fisher



Matthew Rhyner, PhD
Lena Lee
Ara Kulhanjian

Thank You!

Contact us:

Flownamics, Inc.

3025 South Stoughton Road
Madison, WI 53716

phone: 800.932.6989
608.240.1604

e-mail: contact@flownamics.com

web: www.flownamics.com

